



Smart products
for better therapies



**SWISS DENTAL MANUFACTURE WITH GLOBAL EXPERTISE
FOR MORE THAN 75 YEARS.**

At PD, we are focused on providing dental professionals with better and safer solutions to simplify dental therapies. Our innovative dental solutions are a blend of smart yet simple products aimed to improve dental care and revolutionize the industry as we know it.

We are dedicated to promote PRE-ENDO, ENDO and POST-ENDO solutions to help professionals keep the natural teeth of their patients much longer.

Because with the
right tools, anything
is possible.

Yann Gehrig, CEO



History

Founded in 1940 by William Gehrig, Produits Dentaires SA has advanced through the decades, thanks to reputed know-how, outstanding quality and innovative mindset. Today, the company is still family-owned and the third generation strives daily to implement actions answering the challenges of the international dental industry.

Products

Through a comprehensive and smart product range, Produits Dentaires SA offers predictable, safe and efficient dental solutions to improve and simplify PRE-ENDO, ENDO and POST-ENDO dental therapies.

Distribution

Present in more than 100 countries on 5 continents, Produits Dentaires SA has created, over the years, an extensive distribution network of specialized dental dealers and depots. Thanks to our partners, "PD" brand and products have gradually gained international repute, representing Swiss Quality at its best.

Social Responsibility


As a human-oriented and independent family business, Produits Dentaires SA has always considered staff integration a priority in its corporate culture. We also actively support local cultural events and international humanitarian associations through donations and free goods.

Ecological Commitment


Produits Dentaires SA is more than ever aware of the environmental issues related to the production and distribution of its products throughout the world. Compliance with requirements and constant optimization of its processes show the willingness of the company to fulfill its duties towards our environment.


Innovation

Produits Dentaires SA invests in the development of new products which are based on the latest knowledge. As a Swiss manufacturer, we strive to introduce true innovations in line with our high-quality standards in order to facilitate the daily work of dental professionals. In addition, existing products are continuously screened and improved. To meet this goal, Produits Dentaires SA works closely together with national and international expert teams consisting of dentists, dental hygienists and other specialists from the medical field. In addition, several research projects are running in close cooperation with universities and dental schools in Switzerland and worldwide.



Smart therapies to keep natural teeth for a lifetime





9 . PRE-ENDO THERAPY

8 . PROPHYLAXIS

Washing
 Disclosing
 Cleaning
 Polishing
 Home care

23 . ENDO THERAPY

22 . ENDODONTICS

Access
 Desobturation
 Lubrication
 Cleaning
 Irrigation
 Drying
 Medicamentation
 Temporary sealing
 Temporary filling
 Dentin repair
 Final obturation

53 . POST-ENDO THERAPY

52 . RESTORATION

Impression
 Separation
 Contouring
 Final restoration
 Control

69 . SPECIAL SUPPLIES

68 . CLINIC

Prophy paste carrier
 Dispensers
 Cotton
 Patient protection
 Suction
 Cleaning
 Identification

78 . LAB

Soldering

81 . REFERENCES

. INDEX



Prophylaxis

Dental prophylaxis is one of the most widely administered procedures in the dental office. A dental prophylaxis treatment consists of cleaning and polishing the teeth with a paste applied with a prophylactic cup. The ultimate goal of dental prophylaxis is the removal of stain from the teeth. Dental prophylaxis also serves an important role in the prevention of gingivitis and periodontal disease. While infrequent prophylaxis treatments alone are not sufficient to provide an anti-caries benefit, dental prophylaxis as part of a professionally-administered oral hygiene program has been shown to reduce the risk of caries.



PRE-ENDO THERAPY

Root canal infections are primarily caused by deep tooth decay. Fortunately, cavities but also gums diseases can be avoided by proper daily dental hygiene and regular professional dental cleaning. PRE-ENDO therapy by PD consists in advanced and refined dental products specifically designed for efficient and cost-effective preventive dentistry in the dental office and at home. Preserving healthy teeth and gums through, ideally, the most minimally invasive dental treatments is the goal. For PD, it's not about how much dental professionals can get done in a day, but how positively they can impact the dental health of their patients. We are dedicated to developing products that are as simple and efficient as possible, addressing the problem while successfully preserving healthy teeth and gums.

*I have chosen
for me and my dental
prophylaxis team the single use
PD™ ProphyCup and PD™ ProphyPaste.
Both are products I recommend for their
excellent workability, efficiency and cost
effectiveness. I'm very pleased to have
highly satisfied dental hygienists
and patients.*

*Dr Philippe Métillon
Switzerland*

PD™ Effervescent mouth wash tablets

Powerful mouth freshener

Effective effervescent cleansing and refreshing action mouthwash (gargle) tablets. The tablets contain Sodium Bicarbonate and added surfactants all effective for removal of loose dental tissues and debris during and after dental treatment. In addition, the enhanced formulation is ideal for prophylactic home treatment.

Features and benefits

- Peppermint flavor for intense fresh taste
- Excellent value – 1 tablet per glass of mouthrinse is enough
- Economic packing makes them ideal for today's dental practices and home use

Jar of 250 up to 1,000 tablets.

See references on page 81.



PD™ Disclosing tablets

Reveals harmful dental biofilm

Simple solution to educate patients and increase their compliance while simplifying professional tooth cleaning and treatment control. The red stains highlight harmful biofilm and show where cleaning is needed. Furthermore, the use of disclosing tablets is ideal for home care to periodically monitor dental biofilm.

Features and benefits

- Harmless red food dye tablets leaving a mild taste of mint
- Fade quickly from teeth simply by rinsing while allowing the professional to identify problem areas in need of deep cleaning and instruct the patient accordingly
- Economic packing makes them ideal for today's dental practices and home use

Jar of 250 up to 1,000 tablets.

See references on page 81.



PD™ ProphyCup Fire . Firm/Red . Webbed

Advanced way to clean

The unique flare and flex are designed to help interproximal access and dental biofilm removal while the firm body is ideal for severe discolorations. The exclusive design allows teeth contouring, contact and helps in selective stain removal while providing comfort of use.

Features and benefits

- Developed to achieve best tooth cleaning and removal of deposits such as biofilm retentions, calculus residues and young calculus
- 100 % latex-free and made of medical grade biocompatible elastomer to offer extremely smooth and flexible edges for atraumatic treatments
- The strategic 7-webbed design is developed to optimize the dosage of the paste
- Designed with a smart wider base and a webbed interior combined with an extra smooth exterior for best treatment results
- Manufactured with a unique latch mandrel made of fiberglass reinforced plastic

Box of 8 or 48 pieces.

See references on page 81.



PD™ RubberCup . Firm/White . Webbed

Cleaning best seller

The 6-webbed design is developed to optimize the dosage of the paste.

Features and benefits

- Designed to flex and flare around the tooth contours and deliver prophylaxis paste to the tooth surface

Box of 12 pieces.

See references on page 81.



PD™ RubberCup . Firm/White . Serrated

Maximum contact for cleaning

Designed to follow the natural contours of the teeth with flex and flare to reach lingual and interproximal areas.

Features and benefits

- Designed for surface contact, low splatter and stain removal

Box of 12 pieces.

See references on page 81.



PD™ ProphyPaste Pro-F^{1,23} . Medium . Fluoride**Perfect starter for discolored surfaces**

Developed to remove thick plaque and stains while leaving a fresh and smooth minty taste and feel after each treatment. Recommended for use in combination with the exclusive PD™ ProphyRing Set (see page 70).

Features and benefits

- Medium-green, stain remover RDA 170 = 70 microns, gluten-free
- Contains 1.23 % fluoride
- Packed in 2 g single dose disposable cups for easy use and improved hygiene
- Special formula developed for comfort during the polishing procedure
- Developed to rinse easily

Box of 125 pieces.

See references on page 81.



PD™ Pumice powder

Economic and ideal for heavy stains

Air-screened peppermint pumice powder to clean natural teeth and metal fillings.

Features and benefits

- Heavy stain remover 6/0 extra fine = 30 microns
- Economic

Jar of 500 g or 1,000 g.
See references on page 81.



PD™ ProphylPaste Pro-N¹⁰⁰. Medium

Fluoride-free universal paste

Developed to remove thick plaque and stains while leaving a fresh and smooth minty taste and feel after each treatment. The paste without fluoride is suggested prior to composite bonding, sealants and orthodontic bracket applications. Recommended for use in combination with the PD™ ProphylRing Set (see page 70).


Features and benefits

- Medium-pink, stain remover RDA 170 = 70 microns
- Special formula developed for comfort during the polishing procedure
- Developed to rinse easily
- Economic and without fluoride for patients who may be allergic to additives

Jar of 50 g in pack of 1 or 3.

See references on page 81.





PD™ ProphylCup Snow . Soft/White . Webbed

Advanced way to polish

The unique flare and flex are designed to help subgingival and interproximal access and dental biofilm removal. In addition, the soft body is adapted for mild to moderate discolorations. The exclusive design allows teeth contouring, contact and helps in selective stain removal while providing comfort of use.

Features and benefits

- Developed to achieve best tooth cleaning and removal of deposits such as biofilm retentions, calculus residues and young calculus
- 100 % latex-free and made of medical grade biocompatible elastomer to offer extremely smooth and flexible edges designed to gently massage soft tissues
- The strategic 7-webbed design is developed to optimize the dosage of the paste
- Designed with a smart wider base and a webbed interior combined with an extra smooth exterior for best treatment results
- Participate in creating a more shiny, smooth and clean teeth look
- Manufactured with a latch mandrel made of fiberglass reinforced plastic

Box of 8 or 48 pieces.

See references on page 81.



PD™ RubberCup . Soft/Red . Webbed

Polishing best seller

The soft 6-webbed design is developed to optimize the dosage of the paste for best polishing results.

Features and benefits

- Designed to flex and flare around the natural tooth contours and deliver prophyl paste to the tooth surface

Box of 12 pieces.

See references on page 81.



PD™ RubberCup . Soft/Red . Serrated

Maximum contact for polishing

Designed to follow the natural contours of the teeth with extraordinary flex and flare to reach lingual and interproximal areas.

Features and benefits

- Designed for surface contact, low splatter and gentle polishing

Box of 12 pieces.

See references on page 81.



PD™ ProphyPaste Pro-F^{1,2,3} . Fine . Fluoride

Polishing excellence

Developed to remove light stains and biofilm while leaving a fresh and smooth minty taste and feel after each treatment. Recommended for use in combination with the exclusive PD™ ProphyRing Set (see page 70).

Features and benefits

- Fine-green, polishing RDA 120 = 20 microns, gluten-free
- Contains 1.23 % fluoride
- Packed in 2 g single dose disposable cups for easy use and improved hygiene
- Special formula developed for comfort during the polishing procedure
- Developed to rinse easily
- Designed to give a mirror-like shine to natural teeth

Box of 125 pieces.

See references on page 81.





PD™ ProplyPaste Pro-N¹⁰⁰ . Fine

Fluoride-free universal paste

Developed to remove light stains and biofilm while leaving a fresh and smooth minty taste and feel after each treatment. The paste without fluoride is suggested prior to composite bonding, sealants and orthodontic bracket applications. Recommended for use in combination with the PD™ ProplyRing Set (see page 70).

Features and benefits

- Fine-pink, polishing RDA 120 = 20 microns
- Special formula developed for comfort during the polishing procedure
- Developed to rinse easily
- Economic and without fluoride for patients who may be allergic to additives

Jar of 50 g in pack of 1 or 3.

See references on page 81.

PD™ Dental floss . Silk

Easy and comfortable flossing

This all natural dental floss, 100 % pure silk coated in beeswax is a strong, smooth, shred-resistant dental floss that reaches between the teeth and under the gumline to remove dental biofilm that other instruments alone cannot reach.

Features and benefits

- Texturized pure silk fibers are gentle on inflamed or sensitive gums and clean biofilm away to reduce gingivitis and bleeding but also prevent interdental cavities

Plastic dispenser of 25 m spool or glass tube of 100 m spool.

See references on page 81.



PD™ Dental floss . Nylon

Fine waxed floss for tight spaces

Thin, shred-resistant, lightly waxed dental floss particularly effective in very tight contacts between teeth and under the gumline. Specially designed and manufactured to give durability and spreadability to disrupt biofilm deposits and carry them away.

Features and benefits

- Designed to reduce shredding and remove plaque between-teeth and below the gumline while creating a satin-like finish that is extremely soft and gentle

Plastic dispenser of 25 m spool or glass tube of 100 m spool.

See references on page 81.



PD™ Dental sticks

Sustainable Swiss linden wood

The triangular shape and the micro-thin form of the dental sticks are designed to perfectly fit in the interdental spaces and along the gum line while protecting gums and tooth enamel from being damaged during the cleaning process of removing biofilm and food debris.

Features and benefits

- Single ended dental sticks massage carefully the gums and help prevent problems while the intense peppermint flavor leaves the mouth with a clean and fresh feeling after use
- Designed with an extra-short stub to reduce splice formation

Pack of 5 x 25 sticks.

See references on page 81.





Endodontics

Endodontic therapy, also known as endodontic treatment or root canal therapy, is a treatment sequence for the infected pulp of a tooth which results in the elimination of infection and the protection of the decontaminated tooth from future microbial invasion. Root canals, and their associated pulp chamber, are the physical hollows within a tooth that are naturally inhabited by nerve tissue, blood vessels and other cellular entities. Together, these items constitute the dental pulp. Endodontic therapy involves the removal of these structures, the subsequent shaping, cleaning, and decontamination of the hollows with small files and irrigating solutions, and the obturation (filling) of the decontaminated canals. Filling of the cleaned and decontaminated canals is done with an inert filling such as gutta-percha and typically a eugenol-based cement. Endodontics includes both primary and secondary endodontic treatments as well as periradicular surgery which is generally used for teeth that still have potential for salvage. The general objective of endodontic primary treatment is to retain the treated tooth in a healthy state to act as support for a suitable restoration that will allow the tooth to return to normal function, and the specific goal is to prevent or heal apical periodontitis.





ENDO THERAPY

A root canal treatment is one of the most common dental procedures. This treatment can save natural teeth and prevent the need for dental implants or bridges. ENDO therapy by PD consists in an advanced and complete product range specifically designed for efficient and cost effective endodontic treatments. Because every tooth deserves a second (or third) chance, PD is committed to developing dental solutions that will increase endodontic treatment success rates. Collaboration with leading universities and opinion leaders is an important success factor, which is leading us to groundbreaking new developments.

*The MAP System and
PD™ MTA White are my choice
for managing in a predictable and
conservative way treatments such as pulp
capping, perforation repair, apical plug
procedures and micro-surgery.*

*Dr Riccardo Tonini
Italy*



PD™ Gates-Glidden drills

Details make perfection

PD™ Gates-Glidden drills of 32 mm length can be used in any right-angle handpiece to enlarge the coronal portion of the root canal for bulk removal of tooth structure in an easy and reliable manner. Mostly indicated for crown down technique, step back technique, to remove the lingual shoulder during access cavity preparation and to enlarge the canal orifice.

Features and benefits

- State of the art manufacturing process and crafted of corrosion resistant superior stainless steel to offer outstanding cutting efficiency
- Flame shaped instruments notched for easy size identification
- Available in refill or assortment rack

Rack of 6 pieces from size RA 1 to RA 6.

See references on page 82.

PD™ Peeso reamers

Quality that you know and trust

PD™ Peeso reamers of 32 mm length are high quality endodontic tools which are designed to enlarge the cervical portion of the canal, prepare the canal entrance, and remove gutta-percha for post space.

Features and benefits

- State of the art manufacturing process and crafted of corrosion resistant superior stainless steel to offer outstanding cutting efficiency
- Notched for easy size identification, they can be used in any right-angle hand piece
- Sharp edges
- Available in refill or assortment drawer boxes

Rack of 6 pieces from size RA 1 to RA 6.

See references on page 82.



PD™ Eugenate desobturator

Solvent for endodontic filling materials

PD™ Eugenate desobturator is a softening solution made of orange oil that helps remove eugenate base cements and gutta-percha. The softening power is an alternative to power instruments during desobturation.

Features and benefits

- Softening with one drop
- One minute activation
- Natural orange flavoring for patient and user comfort

Bottles of 15 ml.

See references on page 82.



PD™ Resin remover

Solvent for endodontic resin-based fillings

The PD™ Resin remover is a softening solution made of citronella oil that helps remove resin-based root canal fillings. The softening power is an alternative to power instruments during desobturation.

Features and benefits

- Softening with one drop
- One minute activation
- Natural citronella flavoring for patient and user comfort

Bottles of 15 ml.

See references on page 82.



PD™ EDTA 17 % . Gel

Efficient, predictable and safe outcomes

Intended for chemo-mechanical root canal preparation, EDTA helps to remove calcium salts and pulp from decalcified root canals. It facilitates the enlargement of the root canal during instrumentation by softening the walls of the canal. Can be applied as often as necessary during instrumentation and can be removed simply by irrigating with sodium hypochlorite.

Features and benefits

- Gel helps lubricate endodontic instrument and makes penetration easier
- 17 % EDTA helps in chelating the calcium salts and removing smear layer
- Disposable syringe tips to eliminate cross-contamination

Syringe of 3 g.

See references on page 82.



PD™ Root canal barbed broaches

Easy and safe pulp tissue removal


PD metal-handle Barbed Broaches are designed to remove efficiently pulp tissue during root canal procedures.


Features and benefits

- Broaches are made of superior stainless steel
- Flexible barbs are arranged spirally on the cone shaped core for bulk removal of pulp tissue
- Barbs are resilient and do not break in the canal if used properly
- 21 mm long and color coded to ensure proper ISO size selection
- Available in refill or assortment rack

Rack of 12 pieces from size ISO 25 to ISO 60.

See references on page 82.





IrriFlex®

The new standard in dental irrigation

IrriFlex® is a unique irrigation needle made of plastic with a back-to-back 2-side vent design delivering solutions as close as possible to the apex. IrriFlex® is easy to use and delivers fluids in the entire canal system while being a closed-ended needle. Thanks to a soft polypropylene body, IrriFlex® offers an ideal flexibility and reachability for improved irrigation therapies.

Features and benefits

- 30G, 4 % taper, 27 mm long, depth marks (22 mm, 20 mm, 19 mm, 18 mm), luer lock connection
- Soft polypropylene body to make the needle curve and flex easily in complex root canal anatomies
- Smooth glide to allow comfortable solution activation with back and forth movements
- Back-to-back 2-side vent design to secure a balanced irrigation volume and a better control on the root canal irrigation
- Needle fully transparent to improve user comfort and liquid delivery precision
- 60° angulation (cuff-needle) to improve accessibility and comfort
- Plastic body to make the needle atraumatic for hard and soft tissues
- Slightly pre-bendable to meet specific irrigation requirements

Boxes of 20 pcs and 40 pcs.

See references on page 82.



PD™ Plastic pipettes . 1 ml

Precise liquid delivery in narrow spaces


PD ultrathin disposable dental pipettes are used to safely and precisely dispense solutions during endodontic procedures at a controlled rate and reduce the need for cumbersome syringe and needle use.

Features and benefits

- 1 ml capacity, ultrathin tip, color coded and flexible
- Excellent value alternative in terms of quality and price

Box of 5 pieces.

See references on page 82.



PD™ Plastic pipettes . 0.3 ml

Smart and cost-effective liquid delivery

Ultrathin dental pipettes from PD are disposable and a great way to avoid cross-contamination and ideal to dispense low volume solutions safely during endodontic and restorative therapies. These pipettes reduce the need for cumbersome syringe and needle use.

Features and benefits

- 0.3 ml capacity, color coded and flexible
- Excellent value alternative in terms of quality and price

Box of 5 pieces.

See references on page 82.



PD™ EDTA 17 % . Solution

Eases instrumentation and smear layer removal

Intended for chemo-mechanical root canal preparation, EDTA 17 % solution is an effective calcium binding or chelating agent used to facilitate instrumentation of root canals and for smear layer removal. The solution decalcifies the canal walls making it easier to enlarge and shape the canal with files and reamers.

Features and benefits

- Softens the tissues and makes the mechanical drainage of the root canal easier
- Effective chelating agent leaving the dentinal tubules clear of inorganic debris

Bottle of 15 ml.

See references on page 82.





PD™ Paper points . 2 %

Safe and super-absorbent

Sterile and designed in a 2 % (0.02) taper with a 28 mm length to fit precisely prepared root canals and offer a quick drying and efficient moisture removal. PD™ Paper points are binding agent free, highly absorbent, rigid and flexible. Accurate depth measurement and insertion with length marks at 16 mm, 18 mm, 19 mm, 20 mm, 22 mm and 24 mm.

Features and benefits

- Uniformly rolled to extremely close tolerance
- Guaranteed consistency in size and shape
- Highly absorbent and sterilized by gamma-ray
- Accurate insertion of the working length
- White tip to monitor cleaning and drying
- Top color coded to prevent mixing up of sizes
- Long shelf life and freshness assured
- Available in refill or assortment cellpack and drawer box

Cellpack or drawer box of 200 pieces from size ISO 10 to ISO 140.
See references on page 82.

PD™ Paper points . 4 % and 6 %

Greater taper ISO sizes

Sterile and designed in a 4 % (0.04) and 6 % (0.06) greater taper with a 28 mm length to fit precisely prepared root canals and offer a quick drying and efficient moisture removal. PD™ Paper points are binding agent free, highly absorbent, rigid and flexible. Accurate depth measurement and insertion with length marks at 16 mm, 18 mm, 19 mm, 20 mm, 22 mm and 24 mm.

Features and benefits

- Uniformly rolled to extremely close tolerance
- Guaranteed consistency in size and shape
- Highly absorbent and sterilized by gamma-ray
- Accurate insertion of the working length
- White tip to monitor cleaning and drying
- Top color coded to prevent mixing up of sizes
- Long shelf life and freshness assured
- Sold in easy-to-dispense drawer box in counts of 100 paper points
- Available in refill or assortment drawer boxes

Drawer box of 100 pieces from size ISO 15 to ISO 60.
See references on page 82.



PD™ Paper points . Conventional

Ideal for more conically-shaped root canals

PD™ paper points (conventional) have a length of 28 mm and a wider taper design to reduce the number of paper points required to dry a root canal when compared with ISO sizes. Manufactured without any adhesives, paper points are highly absorbent, rigid and flexible.

Features and benefits

- Ideal for drying more conical root canals
- Uniformly rolled to extremely close tolerance
- Guaranteed consistency in size and shape
- Highly absorbent and sterilized by gamma-ray
- Long shelf life and freshness assured
- Sold in easy-to-dispense drawer box in counts of 200 paper points
- Available in refill or assortment drawer boxes

Drawer box of 200 pieces from size XXF to XL.
See references on page 83.





OPACAL

Pre-mixed Calcium Hydroxide paste

Opacal is a creamy, radiopaque, ready-to-use calcium hydroxide paste with a pH above 12.5 recommended for root canal treatments and periapical therapies. The aqueous non-setting formula is ideal for temporary root canal dressings. In addition, the unique paste consistency combined with the special ergonomic syringe and tip design enables to place the paste as close as possible to the apex.

Features and benefits

- High pH value
- High radiopacity
- Ideal creamy consistency
- Perfect and comfortable paste delivery
- Easy to remove (non-setting material)
- Aqueous composition
- Excellent value for money

Syringe of 3 ml with 8 applications tips and 1 finger grip.

See references on page 83.



Pro-Pul-Pan

Reinforced Zinc Oxide Eugenol (ZOE) Cement

Pro-Pul-Pan is an intermediate ivory shade restorative material for class I and class II restorations which is composed of zinc oxide and eugenol reinforced with polymers for temporary filling and sealing. Pro-Pul-Pan may also be used as a base under cements and restorative materials that do not contain resin components, such as amalgams, inlays and onlays.

Features and benefits

- Radiopaque
- Easy to remove
- Simple handling and reliable quality

Box includes 30 g powder and 15 ml liquid.

See references on page 83.



PD™ Zinc Oxide . >99 % pure

Extra pure zinc oxide powder

PD™ Zinc Oxide Powder is designed to formulate temporary filling and sealing materials. Mixed with Eugenol PD™ Zinc Oxide can be used for indirect pulp capping. Furthermore, ZOE bases are compatible with hard and soft tissue.

Features and benefits

- Extra fine powder
- No fillers or additives
- Makes a smooth paste with eugenol
- Provides sufficient compressive strength
- Forms a dressing or lining under non-composite resin

Jar of 100 g, 250 g, 500 g and 1,000 g.
 See references on page 83.



PD™ Eugenol . >99 % pure

Extra pure clove oil

Eugenol oil to be used as a dental dressing with Zinc Oxide. Also indicated for temporary obturation and indirect pulp capping when mixed with Zinc Oxide.

Features and benefits

- Chemically pure with increased effectiveness
- Forms a uniform paste consistency when mixed with Zinc Oxide
- Mixed with Zinc Oxide forms a dressing or lining under cements that do not contain resin components

Bottle of 15 ml, 30 ml and 100 ml.
 See references on page 83.





TEMPOTEC

Eugenol free temporary filling material

TEMPOTEC is a radiopaque self-cured temporary sealing paste used for temporary filling of cavities. It is simple to apply with a filling instrument and provides quick curing in a moist environment. TEMPOTEC is recommended for temporary filling after endodontic treatments.

Features and benefits

- Ready to use
- Radiopaque

Jar of 40 g.

See references on page 83.

MAP One

Endodontic obturations and repairs made simple

The MAP System provides a unique and efficient method for precisely placing MTA by orthograde or retrograde approaches. This professional system is smart and considered the product of choice for perforation, root-end filling, pulp-capping and retro-obturation therapies. The NiTi Memory Shape Needle N° 1 (Ø 1.10 mm) available with MAP One can be manually shaped and adjusts easily to the shape of the root canal for effective placement of mineral trioxide aggregate repair material. The NiTi Memory Shape Needle is by design universal and ideal for both orthograde and retrograde therapies. In combination with PD™ MTA White, the exclusive mineral trioxide aggregate cement developed to be placed with the MAP System, dental professionals have the perfect, simplest and most durable solution for successful and predictable endodontic therapies.

Features and benefits

- Solution for simple, precise and predictable MTA placement
- Extremely high-quality level of material and manufacturing process
- Cost effective method thanks to a faster placement and reduced volume of material needed

Box with 1 syringe, 1 needle, 8 plungers and 1 curette.

See references on page 83.





MAP Needles

Precision for all indications

NiTi Memory Shape needles are by design universal and ideal for both orthograde and retrograde therapies while Classic curved needles are designed for very precise non-surgical endodontic procedures. The Surgery triple-angle needles are best indicated for surgical endodontic procedures and offer excellent access and visibility to the operating field. The needles are manufactured with corrosion resistant alloy. Classic and Surgery needles are available in 0.90 mm and 1.10 mm diameter while NiTi Memory Shape needles are available in 0.90 mm, 1.10 mm and 1.30 mm.

Features and benefits

- Wide range to cover all clinical indications
- Color coded for simplified diameter identification
- Extremely high-quality level of material and manufacturing process

Box with 1 needle, 8 plungers and 1 curette.

See references on page 83.

MAP Intro Kit

Complete solution when precision is a must


MAP Intro Kit is the optimal solution for all professionals looking for flexibility and precision. The memory shape needles available with the kit can be manually shaped to provide perfect access during orthograde and retrograde procedures. The needles take back their initial straight shape during sterilization.

Features and benefits

- Complete solution for simple, precise and predictable MTA placement
- Extremely high-quality level of material and manufacturing process
- Cost effective method thanks to a faster placement and reduced volume of material needed

Box with 1 syringe, 3 needles, 24 plungers, 3 curettes, 1 well, 150 O-rings.

See references on page 83.



MAP Universal Kit


Orthograde and retrograde solution for frequent users

MAP Universal Kit is the complete solution for frequent users of the MAP System and those who need maximum access and placement options. The curved needles available with the kit enable an extremely precise gesture and placement of the material via orthograde access.

Features and benefits

- Complete solution with Classic needles for improved haptic feedback
- Extremely high-quality level of material and manufacturing process
- Cost effective method thanks to a faster placement and reduced volume of material needed

Box with 1 syringe, 5 needles, 40 plungers, 3 curettes, 1 well, 150 O-rings.
See references on page 83.



MAP Surgical Kit

The ultimate retrograde solution for specialists

MAP Surgical Kit is a specific option for the surgery specialists managing complex retrograde obturation cases. The triple-angle needles available with the kit enable easy access and perfect vision of the entire surgical field during complex retrograde obturation procedures.

Features and benefits

- Solution for retrograde therapy specialists
- Extremely high-quality level of material and manufacturing process
- Cost effective method thanks to a faster placement and reduced volume of material needed

Box with 1 syringe, 6 needles, 48 plungers, 2 curettes, 1 well, 100 O-rings.
See references on page 84.





MAP Accessories

Differences are in details

Plastic Plungers are single use and color coded for easy identification. The Paste Cleaning O-rings are color coded for easy identification and delivered with a customized key for an ergonomic and easy removal of excess material around the needle. Finally, the Cleaning Curettes are sterilizable and made of NiTi with a sharp tip to help remove residues of material inside the needles.

Features and benefits

- All accessories are designed to increase durability of needles and improve treatment procedures
- Accessories are color coded for easy identification and improved workflow

16 plungers, 50 O-rings and 1 curette available separately.

See references on page 84.







PD™ MTA White

Biocompatible obturation and repair material

PD™ MTA White is a mineral trioxide aggregate repair cement formulated with calcium tungstate and designed to mix quickly, smoothly, and easily into the right consistency for a variety of procedures. Made of fine hydrophilic particles that set in the presence of water, PD™ MTA White seals off pathways between the root canal system and surrounding tissues, significantly reducing bacterial migration. Its excellent compatibility with the dentinal wall creates a perfect seal that promotes a normal healing response within the root to secure successful procedures and reduce the number of retreatment procedures necessary. PD™ MTA White is indicated for root-end fillings, repair of root canals as an apical plug during apexification, repair of root perforations during root canal therapy, as a consequence of internal resorption, as a pulp capping material and for pulpotomy of primary teeth in the child and adolescent pediatric patient populations. In combination with the MAP System dental professionals secure a precise placement for effective and predictable treatment results.

Features and benefits

- Bismuth oxide free to avoid staining and discoloration
- Optimized particle size allowing homogeneous and complete wetting during mixing
- Avoids bacterial migration and penetration of tissue fluids into the root canal
- Simple manipulation allows to tailor its consistency to suit clinical needs
- Resistant to washout
- Radiopaque
- Perfect marginal sealing
- Non resorbable
- Smart individual dose for optimal use

Box of 4 x 140 mg powder and 3 ml pure H₂O. *See references on page 84.*

EssenSeal®

Ready to become the new benchmark

EssenSeal® is a new generation zinc oxide eugenol sealer featuring tea tree essential oil (Melaleuca) and composed of a powder base and a liquid catalyst. Besides the new pleasant refreshing scent, the unique formulation improves the sealer consistency, the three-dimensional sealing capabilities and the comfort of use. Designed to simplify obturation procedures EssenSeal® mixing ratio is as simple as 1 drop of liquid for 1 spoon of powder. In addition, EssenSeal® is highly radiopaque and non-resorbable. In combination with PD™ Gutta-percha points dental professionals are provided with a perfect and predictable solution for root canal obturations.

Features and benefits

- Pleasant scent
- Convenient powder to liquid ratio
- Mixes easily
- Ideal creamy consistency to facilitate placement
- Small particle size to penetrate easily root anatomy and dentin micro structure
- Three-dimensional sealing
- Radiopaque
- White cement for good aesthetic results
- Highly stable and durable
- Non resorbable
- Can be removed easily
- Working time 60 minutes
- Setting time < 3 hours

Box includes 15 g powder and 10 ml liquid.
See references on page 84.





PD™ Gutta-percha points . 2 %

High quality for optimal results

Gutta-percha is the time-proven fill for root canal therapies. The advanced PD™ Gutta-percha points are designed in a 2 % (0.02) taper and 28 mm length to fit precisely prepared root canals and can be used with either cold obturation techniques such as single cone and lateral compaction or as a master cone for warm vertical compaction. The gutta-percha cones are machine-rolled with precision from the highest quality gutta-percha to ensure true size and optimal treatment results. PD™ Gutta-percha points are flexible yet stiff enough to provide ease during insertion and ideal workability. In addition, the top color code prevents mixing up of sizes and length marks at 16 mm, 18 mm, 19 mm, 20 mm, 22 mm and 24 mm secure an accurate depth measurement and insertion for successful treatments.

Features and benefits

- Uniformly rolled from the finest gutta-percha to extremely close tolerance
- Highly radiopaque
- Accurate insertion of the working length
- Time saving with safe, easy and fast filling
- Top color coded to prevent mixing up of sizes
- Long shelf life and freshness assured
- Available in refill or assortment drawer boxes

Drawer box of 120 pieces from size ISO 10 to ISO 140.

See references on page 84.

PD™ Gutta-percha points . 4 % and 6 %

Greater taper ISO sizes

The advanced PD™ Gutta-percha points with greater taper can be used with either cold obturation techniques such as single cone and lateral compaction or as a master cone for warm vertical compaction. The gutta-percha cones are designed in a 4 % (0.04) and 6 % (0.06) taper and 28 mm length to fit precisely prepared root canals and offer a quicker obturation option for root canals with a greater taper. The gutta percha cones are machine-rolled with precision from the highest quality gutta-percha to ensure true size and optimal treatment results. PD™ Gutta-percha points are flexible yet stiff enough to provide ease during insertion and ideal workability. In addition, the top color code prevents mixing up of sizes and length marks at 16 mm, 18 mm, 19 mm, 20 mm, 22 mm and 24 mm secure an accurate depth measurement and insertion for successful treatments.

Features and benefits

- Uniformly rolled from the finest gutta-percha to extremely close tolerance
- Highly radiopaque
- Accurate insertion of the working length
- Time saving with safe, easy and fast filling
- Top color coded to prevent mixing up of sizes
- Long shelf life and freshness assured
- Available in refill or assortment drawer boxes

Drawer box of 60 pieces from size ISO 15 to ISO 60.
See references on page 84.



PD™ Gutta-percha points . Conventional

Ideal for more conically-shaped root canals

PD™ Gutta-percha points (conventional) have a wider taper design to reduce the number of gutta-percha points required to obturate a root canal when compared with ISO sizes. The advanced PD™ Gutta-percha points can be used with either cold obturation techniques such as single cone and lateral compaction or as a master cone for warm vertical compaction.

Features and benefits

- Ideal for obturating more conical root canals
- 28 mm length, pink, XXF-fine – XL-large sizes
- Highly radiopaque
- Long shelf life and freshness assured
- Sold in easy-to-dispense drawer box in counts of 100 gutta-percha points
- Available in refill or assortment drawer boxes

Drawer box of 100 pieces from size XXF to XL.
See references on page 84.



PD™ Acrylic points . 2 %


Condensation-free obturation

Radiopaque PD™ Acrylic points are 29 mm length, stiff but flexible and mechanically produced to meet ISO tolerances. The acrylic points are designed not to break or twist even in the narrowest and most curved canals and ensure a perfect obturation at the apical level. The acrylic points are equally rigid as silver points but like gutta-percha points, they can be removed by a bur at any time from the root canal. Last but not least, they are available color coded or fully pink in different ISO sizes.

Features and benefits

- Time-saving treatment since they do not need any condensation, though granting a tight apical obturation
- Available in refill or assortment packs

Box of 6 tubes of 10 pieces each from size ISO 10 to ISO 140.
See references on page 84.



PD™ Rosen root canal cement

Rickert's paste silver-free

Rosen root canal cement from PD is a zinc oxide eugenol sealer with thymol iodide that has a powder and liquid catalyst. The sealer comes in an ideal powder liquid presentation to adjust the cement viscosity to the obturation technique. Also, PD Rosen root canal cement is a radiopaque and non-resorbable sealing paste for the permanent obturation of the root canal with systematic insertion of points. In combination with PD™ Gutta-percha points and PD™ Acrylic points dental professionals get a complete obturation solution.


Features and benefits

- Silver-free composition to avoid teeth staining
- Dimensionally stable and non-resorbable for durable obturation

Box includes 15 g powder and 15 ml liquid.

See references on page 83.





Restoration

A dental restoration or dental filling is a treatment to restore the function, comfort, integrity, and morphology of missing tooth structure resulting from caries or external trauma. Although the materials and methods of endodontic and restorative treatment modalities have changed, the ultimate goals remained the same. There is a relationship between endodontic treatment and restorative dentistry.

POST-ENDO THERAPY

Placement of a definitive restoration that seals the tooth coronally is vital for the long-run success of root canal treatment. Lack of a seal totally tight to fluids can end in microorganism contamination of the obturated root canal system and failure of the case. In addition to sealing the tooth coronally, the final restoration should reinforce the remaining tooth structure to avoid fracture and help the tooth return to full function. POST-ENDO therapy by PD consists in a proven range of high-quality products engineered for predictable and efficient restorative therapies. PD is committed to developing high grade restorative dental solutions that will increase root canal treatment success rates. We are dedicated to maintaining patients natural smile for a lifetime while simplifying and improving dental procedures.

*Our experience with the
FIBRAPOST & DRILL-A-POST system
has been very positive, they have delivered
excellent results and suit our
procedures extremely well.*

*Prof. Dr Marcel Wainwright
Sweden*

PD™ Gutta-percha stopping . Sticks

Predictable soft tissue management

Gutta-percha stopping from PD is a radiopaque temporization material made from refined gutta percha and ideal for gingival retraction prior impressions.

Features and benefits

- Ideal to help locate thermal sensitivity
- Easy to heat
- Available in white, pink and ivory color

Sticks packed in 30 g or 120 g.
See references on page 85.

PD™ Hold

Adhesive for alginates and elastomers

Hold is a high-quality adhesive from PD that provides a strong adhesive bond and consistent retention of impression materials. In addition, it helps retain impression material to the perforated or non-perforated metallic and plastic impression trays to provide that extra confidence.

Features and benefits

- Super-adhesive
- Easy removal
- Available in liquid or spray with or without gas

Spray or glass bottle from 30 ml to 120 ml.
See references on page 85.



PD™ Anatomical Wooden wedges

Expert design for quality restorations

Anatomically carved Wooden wedges from PD are made of certified Maplewood from sustainably managed Swiss forests. The interdental wedges ensure an optimal matrix adaptation along the entire cervical margin during multi-functional procedures. Wedges from PD are available in several colors: blue, green, yellow, pink and sized from 11 mm (XS) to 15 mm (L). Perfect for every procedure wedges from PD hold filling and restorative materials in the matrix and reduce the formation of overhang, making polishing easier and quicker.

Features and benefits

- Wedges will not chip or splinter and are designed with a square end for easy removal
- Rounded tip to avoid trauma easily passes over gingiva and concave side conforms matrix band to the anatomy of the tooth
- Carefully and individually manufactured with a perfect finishing and color code
- Small size (XS, 11 mm) is ideal for extremely small spaces and pediatric dentistry

Box of 100 or 400 pieces. *See references on page 85.*



PD™ Interdental Wooden wedges

Made from sustainably managed Swiss forests

These proven wedges ensure better access and stable holding of matrix systems during critical restorative procedures. They are specially crafted from Birchwood for best tooth separation.

Features and benefits

- Wedges are available natural or color coded in four sizes
- Special design provides firm holding and easy, safe insertion
- Rounded base of the tip provides smooth inserting between teeth
- Shaped sides conform to interdental space to provide optimal and uniformed lateral pressure on the matrix band
- Made from birch wood which makes the wedge soft and flexible

Box of 100 or 400 pieces.

See references on page 85.



PD™ Vaseline**Dental quality**

Vaseline is recommended to lubricate the inner surfaces of rubber dam to improve sliding the sheet over the contours of the teeth to overcome the contact areas more easily, and close tightly around the cervix of the tooth. The vaseline is ideal for preparing the buccal area, isolating matrices and protecting tissues.

Features and benefits

- Healthy pure petrolatum jelly for dental use
- Stable, non-water soluble, topical lubricant, unscented, hypoallergenic, non-comedogenic

Tube of 15 g.

See references on page 85.



PD™ Tofflemire matrix retainers

Swiss precision at its best

Tofflemire matrix retainers from PD are made of high-quality stainless steel with superior craftsmanship for optimal results and precision. The ergonomic design provides maximum comfort and safety while the matrix bands offer a perfect stability during condensation of restoration material.

Features and benefits

- Non-slip grip premium quality handle
- High standard polish to avoid gingival trauma but also increase durability and corrosion resistance
- High degree of precision and flexibility while conducting clinical procedures
- Available with or without an assortment of 12 bands and in three sizes (junior, senior and universal)

Box includes 1 retainer and 12 assorted bands.

See references on page 85.



PD™ Tofflemire matrix bands

Memory-free matrices


Tofflemire matrix bands from PD are made of high-quality stainless steel with superior craftsmanship for optimal results and precision. Pre-contoured, thin 0.05 mm and extra-thin 0.035 mm stainless steel bands are easy to mold and burnish for proximal contact. Last but not least, the matrices stay put when burnished against adjacent teeth.

Features and benefits

- Anatomically preformed
- Mirror finish for easy removal
- Smooth edges for patient comfort
- Available in a variety of sizes and contours

Box of 12 pieces.

See references on page 85.



PD™ Transparent strips

Light permeable

Transparent strips from PD ensure complete polymerization and are pre-cut to ease and speed Class II or III restorations. The strips are made of a super thin strong plastic with a flexibility that ensures that the band can be drawn tightly and snugly around the tooth giving excellent marginal adaptation.

Features and benefits

- Crystal clear plastic ensures complete polymerization
- For use with light cure or self-cure composites
- Economic package of 100 pcs, in 8 mm or 10 mm wide
- Available in straight or curved shape, thin (0.075 mm) or extra-thin (0.05 mm)

Pack of 100 pieces.

See references on page 85.



PD™ Transparent rings

Thermoplastic matrices

Transparent rings from PD provide excellent visibility to the tooth structure improving user comfort and composite polymerization during core-build and fillings. The 9 mm height rings are extremely flexible, resistant and adjustable to ensure precise anatomical fit around the cervical outline. Moreover, the rings are exceptionally smooth and produce a composite surface requiring minimal finishing and polishing.

Features and benefits

- Transparent rings to avoid discoloration when in contact with etching gels, bonds and primers
- Made of acetate cellulose to naturally resist deformation and make placement easier
- High transparency for light cure materials
- Economic assortment of 100 pcs

Box of 12 pieces or 100 pieces assorted.

See references on page 85.



DRILL-A-POST

The precision drill companion

DRILL-A-POST is unique tapered reamer system with the ability to remove obturators while preparing a post chamber that perfectly corresponds to FIBRAPOST LUX and FIBRAPOST PLUS. Moreover, the drills follow state of the art manufacturing processes and are crafted of corrosion resistant superior stainless steel.

Features and benefits

- Safe and easy to use
- Reliable
- Tapered for FIBRAPOST LUX and FIBRAPOST PLUS

Rack of 1 or 6 pieces.

See references on page 86.





FIBRAPOST LUX

Light conducting – for dual-cure adhesive technique

FIBRAPOST LUX is recommended to retain restorative cores and reinforce endodontically treated teeth. The post is made of translucent longitudinal glass fibers bundled in a proprietary transparent resin matrix to provide outstanding flexural strength and help to preserve the natural lifelike esthetics that are desired when all ceramic restorations are prescribed. With a longitudinal modulus of elasticity and a shear strength close to that of dentine, its mechanical properties are ideal for a homogeneous reconstruction which minimizes the risks of root fracture. In addition, the post is designed with a series of serrations cut into the surface, to provide self-retention regardless of the cementation media used. Moreover, the gentle taper of the FIBRAPOST LUX corresponds to the natural anatomy of the tooth and perfectly matches the post space created by the DRILL-A-POST.

Features and benefits

- Excellent light translucency for efficient light cure deep into the root canal
- Chemically inert and fully biocompatible
- Excellent radiopacity for better position and control
- Dentin-like elasticity for gentle root treatment
- Pre-silanization for secured adhesion
- Easily removed and trimmed
- Available in single cell packaging blister to avoid cross-contamination risks

Cellpack of 5 pieces and box of 12 pieces.

See references on page 86.

FIBRAPOST PLUS

White – for self-cure adhesive technique

FIBRAPOST PLUS feature outstanding radiopacity to help ensure that they are instantaneously identified on a radiograph at placement as well as at follow up appointments. Due to the dentin-like flexural modulus the post flexes with the tooth and transmits stress more evenly along the remaining root structure to help minimize the instance of root fractures compared to metal. In addition, the white, opaque color ensures excellent aesthetic results, without undesirable coloring in the cervical area of the tooth. It is recommended to use the DRILL-A-POST companion drills for best match when preparing the root canal.

Features and benefits

- Chemically inert and fully biocompatible
- Excellent radiopacity for better position and control
- Dentin-like elasticity for gentle root treatment
- Pre-silanization for secured adhesion
- Easily removed and trimmed
- Available in single cell packaging blister to avoid cross-contamination risks

Cellpack of 5 pieces and box of 12 pieces.

See references on page 86.



PD™ Cavity Liner . Light cure

Calcium hydroxide base material

Cavity Liner from PD is a light cure, radiopaque cavity liner and base material specially formulated for use with adhesives, composites and conventional restorative materials.

Features and benefits

- Chemically bonds to adhesive primers, composites and other resin-based materials, and it micromechanically bonds to dentin
- Light paste consistency to facilitate application
- Protects the pulp and dentin from acid-containing cements and restoratives
- Virtually insoluble in water and oral fluids

Syringe of 4 g and 5 applicator tips.

See references on page 86.



PD™ Calcium Hydroxide powder

Ca(OH)₂ chemically pure

Calcium hydroxide powder for use in root canal therapy provides insulation in deep cavities and is also recommended for temporary root filling, pulp capping and pulp protection.

Features and benefits

- Protects vital pulp from infections and further decay
- Easily removed from root canal by simple irrigation

Jar of 50 g.

See references on page 86.



PD™ C&B Zinc Phosphate cement

High strength permanent cement

PD™ C&B Zinc Phosphate cement is designed for the final cementation of crown and bridge restorations and is refined to provide enhanced physical properties. The cement chemically bonds to the tooth structure and metal which provides excellent strength and marginal integrity for long term restorations. Moreover, it is ideal for securing metal inlays, onlays, posts and orthodontic brackets.

Features and benefits

- High compressive strength and good film thickness with adequate radiopacity
- Available in hand-mix powder-liquid in cream white color

Box includes 30 g powder and 15 ml liquid.

See references on page 86.



PD™ Ionover Filling

Qualities of an authentic glass ionomer

PD™ Ionover Filling is a glass-ionomer restorative material which offers excellent adhesion to both the enamel and the dentine. As it prevents the formation of marginal leakage, retentions are not strictly necessary though undercuts provide additional hold to the filling.

Features and benefits

- Available in Universal shade
- Highly flowable, hydrophilic and easy to handle
- Durable and long lasting with low acid erosion
- Wear resistant
- Radiopaque
- Cost-effective

Box includes 15 g powder and 10 ml liquid.

See references on page 86.



PD™ Ionover Luting

Simple handling and reliable quality

PD™ Ionover Luting is a powder/acid mix glass ionomer luting cement for permanent cementation of crowns, bridges, prefabricated or cast posts and orthodontic appliances, inlays and onlays. The glass ionomer luting cement adheres chemically to tooth-substance and non-precious alloys and can be used as a base or lining under composite, amalgam and glass ionomer filling material. Furthermore, due to its high translucency, natural-looking results are accomplished swiftly.

Features and Benefits

- Available in white shade
- Radiopaque
- Cost-effective

Box includes 15 g powder and 10 ml liquid.

See references on page 86.



PD™ PELLA transparent crown forms

Save time and improve esthetics

PD™ PELLA transparent crown forms are peel off crowns which help contour restorative material to match natural dentition. Their ability to help form a crown in both pediatric and adult dentistry saves time and simplifies the crown process. These transparent peel off crowns are ideal for use with chemical and light-cured composites and strip off easily. In addition, the peel off crowns from PD are made of strong, transparent laminated plastic. It is recommended to use Jiffy tubes from PD for a safe and easy PELLA technique.

Features and benefits

- Allows dental professionals to use their preferred composite for a superior esthetic and shade control
- Automatically contours restorative material to match natural dentition
- Strips off easily leaving a smooth surface
- Thin interproximal walls allow for tight contacts
- Provides sufficient strength for easy handling
- Available in 39 different sizes (right and left pairs) in either assortment box or refills
- Pella crowns assortment is packaged in a unique tray that provides convenient easy storage

Box of 39, 18 or 6 pairs.

See references on page 87.



PD™ Jiffy tubes

Safe and easy dispensing tips

PD™ Jiffy tubes are available straight and curved to aid in the placement of dental materials in difficult access areas. Jiffy tubes are ideal to use with PELLA transparent crown forms.

Features and benefits

- Jiffy tubes are transparent, resistant and flexible
- Available assorted

Box of 24 curved or 36 straight.
See references on page 87.



PD™ Steel abrasive strips

High quality finishing and polishing


PD™ Steel abrasive strips are very flexible and adapt to the contoured proximal regions for all types of filling. They are ideal when used on composites, acrylics, glass ionomer cements and amalgam (metals), and are available with 2, 3, 4, 6 and 8 mm width, in fine, medium or coarse grit. Moreover, they are made double or single faced to fit every specific procedure.

Features and benefits

- Designed to produce a smooth restoration surface to help prevent interproximal plaque buildup for superior finishing and polishing results
- Precision-coated with an aluminum-oxide grit
- Long-lasting, resistant and flexible

Pack of 12 pieces.

See references on page 86.



PD™ Articulating paper

Quick, easy and economical occlusal marking paper

Articulating paper from PD is a diagnostic tool used to highlight occlusal contacts and the distribution of occlusal forces. The double-sided articulating paper from PD is made from premium material to ensure reliable markings every time and the extremely thin, elastic and tear resistant 50- and 80-microns papers guarantees absolute pin-point precision in marking any high spot, occlusion and articulation without smudging. The articulating papers are available thin or thick in I- (strips), C- (half horseshoe) and U- from (horseshoe) but also in one-tone (blue or red) and two-tone (blue and red) for convenience of use.

Features and benefits

- Super-sensitive, non-stick coating marks clearly and accurately in wet or dry fields
- Allows simultaneous recording of the upper and lower arches
- Allows occlusion double check

Box of 12 books of 12 sheets each or 6 books of 12 sheets each.

See references on page 87.



Clinic

Workflow and reliability are important issues in the dental practice. To this end, PD offers a range of proven as well as exceptional accessories and solutions to improve dental procedures, comfort and maintenance.

Lab

The soldering technique is an important aspect of dental technology work. Soldering is the most common technique used to join metal materials. In order to create a bond by soldering, the materials involved must be processed with the highest level of care. To this end, PD offers a high-quality universal solder as well as the dedicated top notch accessory to make every soldering procedure a success.



SPECIAL SUPPLIES

The product range of special supplies at PD is held in high esteem worldwide by thousands of dental professionals. For over 75 years, PD has been recognized for Swiss Quality, precision and advanced chemical and mechanical know-how.

*A range of smart products combining
ergonomics and efficiency for
everyday dentistry.*

*Dr Rémy Balthazard
France*



PD™ ProphyRing Set

A must for polishing in style and comfort

PD™ ProphyRing Set is designed to improve user comfort and tooth cleaning workflow. ProphyRings are compatible with PD™ Disposable empty cups, PD single dose ProphyPaste and most single dose prophy pastes in the market. Stylish and extremely comfortable, ProphyRing Set by PD is the perfect prophylaxis companion.

Features and benefits

- Practical and comfortable alternative to dispense ProphyPaste
- Smart and complete solution for an efficient, ergonomic and hygienic prophylaxis
- Sterilizable ProphyRings

Box of 100 disposable cups and 2 ProphyRings.

See references on page 88.

P-One

The innovative must for dental rooms

P-One is an innovative pellet dispenser for every practice that uses cotton pellets for dental procedures. This unique dispenser has a pop-up system that pushes the cotton pellets to the top with ease. Designed with two different screened grids for easy removal of the desired cotton pellet size one pellet at a time. The large grid is used for PD™ Cotton pellets size 1 and 2 while the small grid is used for PD™ Cotton pellets size 3 to 5. The holder itself is white, and has a sliding top that secures the pellets in place. The body of the P-One pellet dispenser has a stainless steel non-skid bottom for better grip on dental unit desktops. The P-One pellet dispenser provides a safer, cleaner environment for the dental assistant and patient.

Features and benefits

- No more unsuccessful grabbing for cotton pellets
- Easy to assemble
- Innovative design
- Gravity-based expulsion system
- High chemical resistance materials
- Simple to use and easy to disinfect
- Dispenses pellet sizes 1 through 5
- Interchangeable size grid

Box includes 1 dispenser, 1 larger grid and 1 small grid.

See references on page 88.





PD™ Dispenser for cotton pellets

The original cotton pellet dispenser

PD cotton pellet dispenser is very stable and made of high-grade stainless steel with a chrome finish. The tiltop pellet dispenser is designed to allow one-handed operation and represents a hygienic and practical solution for storing and dispensing cotton pellets. The scratch and abrasion-proof chrome finish surface prevents bacteria from adhering and is easy to clean and disinfect. The large grid is used for PD™ Cotton pellets size 1 and 2 while the small grid is used for PD™ Cotton pellets size 3 to 5. Available with a pop-up spring design (model B) for bulk cotton pellets or standard (model A) for refill 3 g boxes.

Features and benefits

- Protects cotton pellets against contamination from spray and aerosols
- Spring design allows availability of the cotton pellets for dispensing without requiring additional action
- PD cotton pellet refill 3 g boxes make restocking easy
- Premium quality
- Extremely durable
- Easy to refill
- Dispenses one pellet at a time
- Hygienic: flip down lid to prevent contamination of pellets
- Sterilizable

Box includes 1 dispenser.

See references on page 88.

PD™ EZ Dropper

The solution for all dental solutions

PD™ EZ Dropper is an elegant chair-side dental solution tray. Complete with dropper bottles and dishes, designed to dispense a variety of solutions accurately, conveniently and economically, while reducing the risk of cross contamination. Available in several product configurations with the choice of #1, #2, #3 or #4 tiltop dispenser for cotton pellets.

Features and benefits

- Anti-reflux pump offers drop-by-drop control of the flow, allowing for economical dispensing of solutions
- Color-coded dropper bottles and dishes allow immediate, accurate identification
- Variety of solutions on one tray saves time and is convenient
- Each tray arrives with two dappen dishes of each color

Box includes 1 tray, drop bottles, dappen dishes and accessories as per selection. *See references on page 88.*



PD™ Dappen dish

Professional dental mixing glass

PD™ Dappen dish mixing glass is made of high quality and durable glass. Designed with a thick and heavy base to remain steady and avoid spills it is suitable for working with small amounts of dental solutions, medicaments and fillings. PD™ Dappen dishes are ideal for mixing dental acrylics, amalgam, resins, sealants and other materials.

Features and benefits

- Premium quality
- Highly durable and crack resistant
- Perfect size for minimizing waste and contamination
- Withstand any liquid used in a dental operator
- Sterilizable

Box includes 1 dappen dish.
See references on page 88.





PD™ Cotton pellets

Super soft and extremely absorbent

PD™ Cotton pellets are produced following high manufacturing standards to ensure the highest precision and quality to meet the need of demanding dental professionals when drying surfaces or applying medicaments. The automated production and quality control associated with the best cotton quality ensure constant high absorbency and uniform size and shape. 100 % cotton pellets are designed to provide outstanding quality and absorbency in a variety of uses. Furthermore, the high-quality PD cotton pellet dispensers keep cotton pellets within reach during treatment protecting them hygienically from spray, dust and cross contamination. The cotton pellet refills are designed for use in the standard tiltop cotton pellet dispenser (model A).

Features and benefits

- High quality raw material
- Highly absorbent
- Uniform sizes
- Eliminates cotton waste
- Available in convenient refills or bulk in sizes 1 through 5

Box of 6 refills of 3 g each from size 1 to 5 or in bulk of 18 g or 32 g.
See references on page 88.

PD™ Napkin holders

Colorful, autoclavable and long lasting

PD sterilizable Napkin and bib holders are available in durable nickel-plated polished chrome with strong metal chain or premium quality medical grade soft plastic, both with alligator clips for secure holding and long life.

Features and benefits

- Sterilizable
- Stainless steel alligator clips
- Bright, attractive colors
- 100 % latex-free

Pack of 1 piece.

See references on page 88.

PD™ Saliva ejectors

Medical PVC with fixed weld copper-wire

PD™ Saliva ejectors are designed for maximum patient comfort and maximum aspiration. The specific compact design tip allows better accessibility in dental procedures. Last but not least, the metal wire-reinforced is fitted solidly in the tube, to avoid any chance of it coming out of either end and will guarantee that the desired shape is maintained perfectly.

Features and benefits

- High flexibility devoid kinks to guarantee high filtration capacity
- Basket strainer tip permits maximum fluid aspiration without clogging
- Vacuum suction vents reduce risk of tissue suction
- One fits all standard hose ends
- Easily forms into desired configuration
- Tip design eliminates suction of soft tissue
- Easy to use with non-rusting alloy wire
- Latex free transparent PVC
- Disposable
- Available in economical packs of 100 pcs

Pack of 100 pieces. See references on page 88.





PD™ Bur cleaning brush

Fast and easy cleaning of burs

PD™ Bur cleaning brush is ideal for debridement of carbide burs, diamonds, endodontic instruments and inserts. The brush removes debris from the cutting flutes, therefore prolonging the life span of burs and providing better cutting efficiency.

Features and benefits

- Flat brush
- Adjustable
- Brass bristles

Pack of 1 piece.

See references on page 88.

PD™ Marking rings

Identify, organize and customize instrumentation

PD™ Marking rings are universal color-coded bands made from silicone for identifying and organizing instruments. They are useful for maintaining proper instrument order in cassettes and identifying instruments by procedure. The rings resist all heat sterilization methods and are latex free.


Features and benefits

- Aid in procedure sequence and operator efficiency
- No sticking or cracking in the autoclave
- Latex free
- Small, medium and large sizes are available in 10 colors

Box of 48 pieces.

See references on page 89.





PD™ Universal solder . Silver base

The most trusted solder pre-filled with fluoride flux

PD™ Universal solder is a fine seamless silver base tube pre-filled with fluoride flux. The silver solder/flux sticks save time and effort by placing the flux exactly where it is needed with the solder and they are dead soft for easy handling and cadmium free for safety. In addition, the low differential between melting and flow points provides excellent control. Last but not least, the silver solder sticks with the flux in the center provide the perfect solder/flux ratio for fine soldering.

Features and benefits

- Fluoride flux for cleaning away oxidizing agents
- Creamy smooth formula flows easily into joints
- Low melting point: 630° C
- For all brazing works
- Cadmium free
- Fine soldering
- Time saving

Pack of 12 sticks.

See references on page 89.

PD™ Micro-torch

Professional hand-held butane refillable micro torch

PD™ Micro-torch is a Piezo Electronic Ignition gas burner for soldering, welding, softening wax, melting gold and fusing plastic. This hand torch features a powerful 1300°C/2500°F butane blue torch flame and air flow adjustment dial for superior flame control with pinpoint accuracy. In addition, the main tank body features a high-quality anodized finish design for a superior durability and extra safe grip.

Features and benefits

- Lightweight and compact
- Instant ignition system
- Precise flame adjustment
- Easy refill

- ### Features and benefits
- Lightweight and compact
 - Instant ignition system
 - Precise flame adjustment
 - Easy refill

Box of 1 device.

See references on page 89.



PRE-ENDO THERAPY

| PD™ Effervescent mouth wash tablets | REF |
|-------------------------------------|-------|
| Jar of 250 tablets | 66301 |
| Jar of 500 tablets | 66302 |
| Jar of 1,000 tablets | 66303 |

| PD™ Disclosing tablets | REF |
|--------------------------------|-------|
| Jar of 100 tablets pink | 11611 |
| Jar of 500 tablets pink | 11613 |
| Jar of 1,000 tablets pink | 11612 |

| PD™ PropHyCup Fire . Firm/Red . Webbed | REF |
|--|-------|
| Box of 8 pcs latch-type | 47061 |
| Box of 48 pcs latch-type | 47062 |

| PD™ PropHyCup Snow . Soft/White . Webbed | REF |
|--|-------|
| Box of 8 pcs latch-type | 47051 |
| Box of 48 pcs latch-type | 47052 |

| PD™ RubberCup . Firm/White . Webbed | REF |
|-------------------------------------|-------|
| Box of 12 pcs latch-type | 46602 |

| PD™ RubberCup . Firm/White . Serrated | REF |
|---------------------------------------|-------|
| Box of 12 pcs latch-type | 46612 |

| PD™ RubberCup . Soft/Red . Webbed | REF |
|-----------------------------------|-------|
| Box of 12 pcs latch-type | 46502 |

| PD™ RubberCup . Soft/Red . Serrated | REF |
|-------------------------------------|-------|
| Box of 12 pcs latch-type | 46512 |

| PD™ PropHyPaste Pro-F ^{1,23} . Medium . Fluoride | REF |
|---|-------|
| Box of 125 pcs single dose cups | 47025 |

| PD™ PropHyPaste Pro-N ¹⁰⁰ . Medium | REF |
|---|-------|
| Jar of 50 g | 47015 |
| Jar of 50 g x 3 | 47016 |

| PD™ PropHyPaste Pro-F ^{1,23} . Fine . Fluoride | REF |
|---|-------|
| Box of 125 pcs single dose cups | 47020 |

| PD™ PropHyPaste Pro-N ¹⁰⁰ . Fine | REF |
|---|-------|
| Jar of 50 g | 47010 |
| Jar of 50 g x 3 | 47011 |

| PD™ Pumice powder | REF |
|--------------------------------|-------|
| Jar of 500 g extra-fine | 46403 |
| Jar of 1,000 g extra-fine | 46404 |

| PD™ Dental floss . Silk | REF |
|---------------------------------|-------|
| Plastic dispenser of 25 m spool | 66012 |
| Glass tube of 100 m spool | 66013 |

| PD™ Dental floss . Nylon | REF |
|---------------------------------|-------|
| Plastic dispenser of 25 m spool | 66112 |
| Glass tube of 100 m spool | 66113 |

| PD™ Dental sticks | REF |
|---|-------|
| Pack of 5 x 25 sticks single ended | 65920 |

ENDO THERAPY

| PD™ Gates-Glidden drills | | | | REF |
|--------------------------|---------|---------------|--|--------------|
| N° | Ø | | | |
| RA 1-6 | assort. | rack of 6 pcs | | 20362 |
| RA 1 | 0.50 | rack of 6 pcs | | 20302 |
| RA 2 | 0.70 | rack of 6 pcs | | 20312 |
| RA 3 | 0.90 | rack of 6 pcs | | 20322 |
| RA 4 | 1.10 | rack of 6 pcs | | 20332 |
| RA 5 | 1.30 | rack of 6 pcs | | 20342 |
| RA 6 | 1.50 | rack of 6 pcs | | 20352 |

| PD™ Peeso reamers | | | | REF |
|-------------------|---------|---------------|--|--------------|
| N° | Ø | | | |
| RA 1-6 | assort. | rack of 6 pcs | | 20462 |
| RA 1 | 0.70 | rack of 6 pcs | | 20402 |
| RA 2 | 0.90 | rack of 6 pcs | | 20412 |
| RA 3 | 1.10 | rack of 6 pcs | | 20422 |
| RA 4 | 1.30 | rack of 6 pcs | | 20432 |
| RA 5 | 1.50 | rack of 6 pcs | | 20442 |
| RA 6 | 1.70 | rack of 6 pcs | | 20452 |

| PD™ Eugenate desobturator | | REF |
|---------------------------|--|--------------|
| Bottle of 15 ml | | 11555 |

| PD™ Resin remover | | REF |
|-------------------|--|--------------|
| Bottle of 15 ml | | 11556 |

| PD™ EDTA 17 % . Gel | | REF |
|---------------------|--|--------------|
| Syringe of 3 g | | 11566 |

| PD™ Root canal barbed broaches | | | | REF |
|--------------------------------|----------|----------------|--|--------------|
| 25-60 | assort. | rack of 12 pcs | | 21157 |
| 25 | xxx-fine | rack of 12 pcs | | 21151 |
| 30 | xx-fine | rack of 12 pcs | | 21152 |
| 35 | x-fine | rack of 12 pcs | | 21153 |
| 40 | fine | rack of 12 pcs | | 21154 |
| 50 | medium | rack of 12 pcs | | 21155 |
| 60 | coarse | rack of 12 pcs | | 21156 |

| IrriFlex® | | REF |
|---------------|--|--------------|
| Box of 20 pcs | | 21400 |
| Box of 40 pcs | | 21401 |

| PD™ EDTA 17 % . Solution | | REF |
|--------------------------|--|--------------|
| Bottle of 15 ml | | 11565 |

| PD™ Plastic pipettes . 1 ml | | REF |
|-----------------------------|-------------|--------------|
| Box of 5 pcs | color coded | 11752 |

| PD™ Plastic pipettes . 0.3 ml | | REF |
|-------------------------------|-------------|--------------|
| Box of 5 pcs | color coded | 11750 |

| PD™ Paper points . 2 % | | | REF CELLPACK | REF DRAWER BOX* |
|------------------------|---------|---------|-----------------|-----------------------|
| 15-35 (15-40)* | assort. | 200 pcs | 25931 | 25908 |
| 10 | | 200 pcs | 25924 | 25901 |
| 15 | | 200 pcs | 25925 | 25902 |
| 20 | | 200 pcs | 25926 | 25903 |
| 25 | | 200 pcs | 25927 | 25904 |
| 30 | | 200 pcs | 25928 | 25905 |
| 35 | | 200 pcs | 25929 | 25906 |
| 40 | | 200 pcs | 25930 | 25907 |
| 40-60 (45-80)* | assort. | 200 pcs | 25938 | 25915 |
| 45 | | 200 pcs | 25932 | 25909 |
| 50 | | 200 pcs | 25933 | 25910 |
| 55 | | 200 pcs | 25934 | 25911 |
| 60 | | 200 pcs | 25935 | 25912 |
| 70 | | 200 pcs | 25936 | 25913 |
| 80 | | 200 pcs | 25937 | 25914 |
| 90-140 | assort. | 120 pcs | | 25922 |

| PD™ Paper points . 4 % and 6 % | | REF 4 % | REF 6 % |
|-----------------------------------|-------------------------------|--------------|--------------|
| 15 | drawer box of 100 pcs | 25947 | |
| 20-45 | assort. drawer box of 100 pcs | 25954 | 25969 |
| 20 | drawer box of 100 pcs | 25948 | 25963 |
| 25 | drawer box of 100 pcs | 25949 | 25964 |
| 30 | drawer box of 100 pcs | 25950 | 25965 |
| 35 | drawer box of 100 pcs | 25951 | 25966 |
| 40 | drawer box of 100 pcs | 25952 | 25967 |
| 45 | drawer box of 100 pcs | 25953 | 25968 |
| 50 | drawer box of 100 pcs | 25955 | |
| 60 | drawer box of 100 pcs | 25957 | |

| PD™ Paper points . Conventional | REF |
|--------------------------------------|--------------|
| XXF-XL assort. drawer box of 200 pcs | 25982 |
| XXF drawer box of 200 pcs | 25976 |
| XF drawer box of 200 pcs | 25977 |
| F drawer box of 200 pcs | 25978 |
| M drawer box of 200 pcs | 25979 |
| L drawer box of 200 pcs | 25980 |
| XL drawer box of 200 pcs | 25981 |

| OPACAL | REF |
|-----------------|--------------|
| Syringe of 3 ml | 11020 |

| Pro-Pul-Pan | REF |
|---|--------------|
| Box includes 30 g powder and 15 ml liquid | 10260 |
| Powder bottle of 30 g | 10261 |

| PD™ Zinc Oxide . >99 % pure | REF |
|-----------------------------|--------------|
| Jar of 100 g | 11431 |
| Jar of 250 g | 11432 |
| Jar of 500 g | 11430 |
| Jar of 1,000 g | 11433 |

| PD™ Eugenol . >99 % pure | REF |
|--------------------------|--------------|
| Bottle of 15 ml | 11420 |
| Bottle of 30 ml | 11421 |
| Bottle of 100 ml | 11422 |

| TEMPOTEC | REF |
|-------------|--------------|
| Jar of 40 g | 10204 |

| PD™ Rosen root canal cement | REF |
|---|--------------|
| Box includes 15 g powder and 15 ml liquid | 10245 |

| MAP One | REF |
|--|--------------|
| Content | 20296 |
| 1 NiTi Memory Shape Needle N° 1 8 Plastic Plungers N° 1 1 Cleaning Curette N° 1 1 Stainless Steel Syringe | |

| MAP Needles | REF |
|--------------------------|--------------|
| NiTi Memory Shape | |
| N° 0 yellow 0.90 mm | 20210 |
| N° 1 red 1.10 mm | 20211 |
| N° 2 blue 1.30 mm | 20212 |

| | |
|---------------------|--------------|
| Classic | |
| N° 0 yellow 0.90 mm | 20220 |
| N° 1 red 1.10 mm | 20221 |

| | |
|----------------------------|--------------|
| Surgery | |
| N° 0 yellow 0.90 mm, right | 20225 |
| N° 0 yellow 0.90 mm, left | 20226 |
| N° 1 red 1.10 mm, right | 20227 |
| N° 1 red 1.10 mm, left | 20228 |

| MAP Intro Kit | REF |
|----------------|--------------|
| Content | 20285 |

1 NiTi Memory Shape Needle N° 0
1 NiTi Memory Shape Needle N° 1
1 NiTi Memory Shape Needle N° 2
8 Plastic Plungers N° 0
8 Plastic Plungers N° 1
8 Plastic Plungers N° 2
1 Cleaning Curette N° 0
1 Cleaning Curette N° 1
1 Cleaning Curette N° 2
1 Stainless steel Syringe
1 Paste Cleaning O-rings N° 0
1 Paste Cleaning O-rings N° 1
1 Paste Cleaning O-rings N° 2
1 Stainless steel Well
1 Aluminum case 185 x 95 x 30 mm

| MAP Universal Kit | REF |
|-------------------|--------------|
| Content | 20286 |

1 NiTi Memory Shape Needle N° 0
1 NiTi Memory Shape Needle N° 1
1 NiTi Memory Shape Needle N° 2
1 Classic Needle N° 0
1 Classic Needle N° 1
16 Plastic Plungers N° 0
16 Plastic Plungers N° 1
8 Plastic Plungers N° 2
1 Cleaning Curette N° 0
1 Cleaning Curette N° 1
1 Cleaning Curette N° 2
1 Stainless steel Syringe
1 Paste Cleaning O-rings N° 0
1 Paste Cleaning O-rings N° 1
1 Paste Cleaning O-rings N° 2
1 Stainless steel Well
1 Aluminum case 185 x 95 x 30 mm



| MAP Surgical Kit | REF |
|----------------------------------|--------------|
| Content | 20287 |
| 1 Surgery Needle N° 0 right | |
| 1 Surgery Needle N° 0 left | |
| 1 Surgery Needle N° 1 right | |
| 1 Surgery Needle N° 1 left | |
| 1 NiTi Memory Shape Needle N° 0 | |
| 1 NiTi Memory Shape Needle N° 1 | |
| 24 Plastic Plungers N° 0 | |
| 24 Plastic Plungers N° 1 | |
| 1 Cleaning Curette N° 0 | |
| 1 Cleaning Curette N° 1 | |
| 1 Stainless steel Syringe | |
| 1 Paste Cleaning O-rings N° 0 | |
| 1 Paste Cleaning O-rings N° 1 | |
| 1 Stainless steel Well | |
| 1 Aluminum case 185 x 95 x 30 mm | |

| MAP Accessories | REF |
|---------------------------------|--------------|
| Plastic plungers | |
| N° 0 yellow pack of 16 plungers | 20200 |
| N° 1 red pack of 16 plungers | 20201 |
| N° 2 blue pack of 16 plungers | 20202 |

| | |
|-------------------------------|--------------|
| Paste cleaning O-rings | |
| N° 0 yellow box of 50 O-rings | 20252 |
| N° 1 red box of 50 O-rings | 20253 |
| N° 2 blue box of 50 O-rings | 20254 |

| | |
|--------------------------|--------------|
| Cleaning curettes | |
| N° 0 yellow | 20240 |
| N° 1 red | 20241 |
| N° 2 blue | 20243 |

| PD™ MTA White | REF |
|---|--------------|
| Box includes 4 x 140 mg powder and 3 ml pure H ₂ O | 20190 |

| EssenSeal® | REF |
|---|--------------|
| Box includes 15 g powder and 10 ml liquid | 10246 |

| PD™ Gutta-percha points . 2 % | REF |
|--------------------------------------|--------------|
| 10 drawer box of 120 pcs | 24001 |
| 15-40 assort. drawer box of 120 pcs | 24008 |
| 15 drawer box of 120 pcs | 24002 |
| 20 drawer box of 120 pcs | 24003 |
| 25 drawer box of 120 pcs | 24004 |
| 30 drawer box of 120 pcs | 24005 |
| 35 drawer box of 120 pcs | 24006 |
| 40 drawer box of 120 pcs | 24007 |
| 45-80 assort. drawer box of 120 pcs | 24015 |
| 45 drawer box of 120 pcs | 24009 |
| 50 drawer box of 120 pcs | 24010 |
| 55 drawer box of 120 pcs | 24011 |
| 60 drawer box of 120 pcs | 24012 |
| 70 drawer box of 120 pcs | 24013 |
| 80 drawer box of 120 pcs | 24014 |
| 90-140 assort. drawer box of 120 pcs | 24022 |

| PD™ Gutta-percha points . 4 % and 6 % | REF 4 % | REF 6 % |
|---------------------------------------|--------------|--------------|
| 15 drawer box of 60 pcs | 24024 | |
| 20-45 assort. drawer box of 60 pcs | 24031 | 24046 |
| 20 drawer box of 60 pcs | 24025 | 24040 |
| 25 drawer box of 60 pcs | 24026 | 24041 |
| 30 drawer box of 60 pcs | 24027 | 24042 |
| 35 drawer box of 60 pcs | 24028 | 24043 |
| 40 drawer box of 60 pcs | 24029 | 24044 |
| 45 drawer box of 60 pcs | 24030 | 24045 |
| 50 drawer box of 60 pcs | 24032 | |
| 60 drawer box of 60 pcs | 24034 | |

| PD™ Gutta-percha points . Conventional | REF |
|--|--------------|
| XXF-XL assort. drawer box of 100 pcs | 24059 |
| XXF drawer box of 100 pcs | 24053 |
| XF drawer box of 100 pcs | 24054 |
| F drawer box of 100 pcs | 24055 |
| M drawer box of 100 pcs | 24056 |
| L drawer box of 100 pcs | 24057 |
| XL drawer box of 100 pcs | 24058 |

| PD™ Acrylic points . 2 % | REF PINK | REF COLOR CODED |
|---------------------------------|--------------|-----------------|
| 10 6 tubes x 10 pcs | 26702 | 26902 |
| 15-40 assort. 6 tubes x 10 pcs | 26709 | 26909 |
| 15 6 tubes x 10 pcs | 26703 | 26903 |
| 20 6 tubes x 10 pcs | 26704 | 26904 |
| 25 6 tubes x 10 pcs | 26705 | 26905 |
| 30 6 tubes x 10 pcs | 26706 | 26906 |
| 35 6 tubes x 10 pcs | 26707 | 26907 |
| 40 6 tubes x 10 pcs | 26708 | 26908 |
| 45-80 assort. 6 tubes x 10 pcs | 26717 | 26917 |
| 45 6 tubes x 10 pcs | 26710 | 26910 |
| 50 6 tubes x 10 pcs | 26711 | 26911 |
| 55 6 tubes x 10 pcs | 26713 | 26913 |
| 60 6 tubes x 10 pcs | 26714 | 26914 |
| 70 6 tubes x 10 pcs | 26715 | 26915 |
| 80 6 tubes x 10 pcs | 26716 | 26916 |
| 90-140 assort. 6 tubes x 10 pcs | 26725 | 26925 |
| 90 6 tubes x 10 pcs | 26718 | 26918 |
| 100 6 tubes x 10 pcs | 26719 | 26919 |
| 110 6 tubes x 10 pcs | 26720 | 26920 |
| 120 6 tubes x 10 pcs | 26722 | 26922 |
| 130 6 tubes x 10 pcs | 26723 | 26923 |
| 140 6 tubes x 10 pcs | 26724 | 26924 |

POST-ENDO THERAPY

| PD™ Gutta-percha stopping Sticks | | REF 30 g | REF 120 g |
|-------------------------------------|----------------|-------------|--------------|
| Assort. | pack of sticks | 10220 | 10225 |
| White | pack of sticks | 10223 | 10228 |
| Pink | pack of sticks | 10222 | 10227 |
| Ivory | pack of sticks | 10221 | 10226 |

| PD™ Anatomical Wooden wedges | | REF 100 pcs | REF 400 pcs |
|---------------------------------|-------------|----------------|----------------|
| 11 mm-15 mm | assort. | 21785 | 21786 |
| 11 mm | box, pink | 21781 | 21798 |
| 12 mm | box, blue | 21782 | 21795 |
| 13 mm | box, green | 21783 | 21796 |
| 15 mm | box, yellow | 21784 | 21797 |

| PD™ Interdental Wooden wedges | | REF 100 pcs | REF 400 pcs |
|----------------------------------|---------------------|----------------|----------------|
| Small-Wide | assort. plain | 21725 | 21735 |
| Small | plain | 21721 | 21731 |
| Medium | plain | 21722 | 21732 |
| Long | plain | 21723 | 21733 |
| Wide | plain | 21724 | 21734 |
| Small-Wide | assort. color coded | 21745 | 21755 |
| Small | plain | 21741 | |
| Medium | red | 21742 | |
| Long | yellow | 21743 | |
| Wide | blue | 21744 | |

| PD™ Vaseline | | REF |
|--------------|--|-------|
| Tube of 15 g | | 10651 |

| PD™ Hold | | REF |
|------------------------------|--------|-------|
| Aerosol spray for alginates | 120 ml | 11131 |
| Spray for alginates | 100 ml | 11132 |
| Liquid for alginates | 30 ml | 11130 |
| Aerosol spray for elastomers | 120 ml | 11135 |

| PD™ Tofflemire matrix retainers | | REF |
|---------------------------------|------------------------|-------|
| Junior retainer | 12 assorted thin bands | 41501 |
| Senior retainer | 12 assorted thin bands | 41502 |
| Universal retainer | 12 assorted thin bands | 41503 |
| Junior retainer | | 41504 |
| Senior retainer | | 41505 |
| Universal retainer | | 41506 |

| PD™ Tofflemire matrix bands | | REF X-TIN | REF THIN |
|--------------------------------|-------------------------------|--------------|-------------|
| N° 1-15 assort. | box of 4x#1, 3x#2, 1x#3-15 | 41518 | 41528 |
| N° 1 adults | box of 12 pcs | 41511 | 41521 |
| N° 2 adults | box of 12 pcs | 41512 | 41522 |
| N° 3 adults | box of 12 pcs | 41513 | 41523 |
| N° 11 adults | box of 12 pcs | 41514 | 41524 |
| N° 13 children | box of 12 pcs | 41515 | 41525 |
| N° 14 children | box of 12 pcs | 41516 | 41526 |
| N° 15 children | box of 12 pcs | 41517 | 41527 |

| PD™ Transparent strips | | REF X-TIN | REF THIN |
|------------------------|------------------------|--------------|-------------|
| Straight | pack of 100 pcs, 8 mm | 42011 | 42012 |
| Straight | pack of 100 pcs, 10 mm | 42021 | 42022 |
| Curved | pack of 100 pcs, 8 mm | 42211 | 42212 |
| Curved | pack of 100 pcs, 10 mm | 42221 | 42222 |

| PD™ Transparent rings | | REF |
|----------------------------|----------------|-------|
| N° 1-20 assort. | box of 100 pcs | 10760 |
| N° 1 (int. circum. 15 mm) | box of 12 pcs | 10701 |
| N° 2 (int. circum. 16 mm) | box of 12 pcs | 10702 |
| N° 3 (int. circum. 17 mm) | box of 12 pcs | 10703 |
| N° 4 (int. circum. 18 mm) | box of 12 pcs | 10704 |
| N° 5 (int. circum. 20 mm) | box of 12 pcs | 10706 |
| N° 6 (int. circum. 21 mm) | box of 12 pcs | 10707 |
| N° 7 (int. circum. 22 mm) | box of 12 pcs | 10708 |
| N° 8 (int. circum. 23 mm) | box of 12 pcs | 10709 |
| N° 9 (int. circum. 25 mm) | box of 12 pcs | 10711 |
| N° 10 (int. circum. 26 mm) | box of 12 pcs | 10712 |
| N° 11 (int. circum. 27 mm) | box of 12 pcs | 10713 |
| N° 12 (int. circum. 28 mm) | box of 12 pcs | 10714 |
| N° 13 (int. circum. 30 mm) | box of 12 pcs | 10716 |
| N° 14 (int. circum. 31 mm) | box of 12 pcs | 10717 |
| N° 15 (int. circum. 32 mm) | box of 12 pcs | 10718 |
| N° 16 (int. circum. 33 mm) | box of 12 pcs | 10719 |
| N° 17 (int. circum. 35 mm) | box of 12 pcs | 10721 |
| N° 18 (int. circum. 36 mm) | box of 12 pcs | 10722 |
| N° 19 (int. circum. 37 mm) | box of 12 pcs | 10723 |
| N° 20 (int. circum. 38 mm) | box of 12 pcs | 10724 |

| DRILL-A-POST | | | | REF |
|--------------|---------|---------------|--|--------------|
| N° 1-6 | assort. | rack of 6 pcs | | 22507 |
| N° 1 | | rack of 6 pcs | | 22501 |
| N° 2 | | rack of 6 pcs | | 22502 |
| N° 3 | | rack of 6 pcs | | 22503 |
| N° 4 | | rack of 6 pcs | | 22504 |
| N° 5 | | rack of 6 pcs | | 22505 |
| N° 6 | | rack of 6 pcs | | 22506 |
| N° 1 | | rack of 1 pce | | 22511 |
| N° 2 | | rack of 1 pce | | 22512 |
| N° 3 | | rack of 1 pce | | 22513 |
| N° 4 | | rack of 1 pce | | 22514 |
| N° 5 | | rack of 1 pce | | 22515 |
| N° 6 | | rack of 1 pce | | 22516 |

| PD™ Cavity Liner . Light cure | | REF |
|--------------------------------------|--|--------------|
| Syringe of 4 g and 5 applicator tips | | 11035 |

| PD™ Calcium Hydroxide powder | | REF |
|------------------------------|--|--------------|
| Jar of 50 g | | 11025 |

| PD™ C&B Zinc Phosphate cement | | REF |
|---|--|--------------|
| Box includes 30 g powder and 15 ml liquid | | 10230 |

| FIBRAPOST LUX | | | REF |
|---------------|-------------------|--------------------------|--------------|
| N° 1 | cellpack of 5 pcs | | 22641 |
| N° 2 | cellpack of 5 pcs | | 22642 |
| N° 3 | cellpack of 5 pcs | | 22643 |
| N° 4 | cellpack of 5 pcs | | 22644 |
| N° 1-4 | assort. | box of 24 pcs, 4 reamers | 22589 |
| N° 1-4 | assort. | box of 12 pcs | 22594 |
| N° 1 | box of 12 pcs | | 22590 |
| N° 2 | box of 12 pcs | | 22591 |
| N° 3 | box of 12 pcs | | 22592 |
| N° 4 | box of 12 pcs | | 22593 |

| PD™ Ionover Filling | | REF |
|---|--|--------------|
| Box includes 15 g powder and 10 ml liquid, C1 | | 10390 |

| PD™ Ionover Luting | | REF |
|---|--|--------------|
| Box includes 15 g powder and 10 ml liquid | | 10385 |

| FIBRAPOST PLUS | | | REF |
|----------------|-------------------|--------------------------|--------------|
| N° 1 | cellpack of 5 pcs | | 22631 |
| N° 2 | cellpack of 5 pcs | | 22632 |
| N° 3 | cellpack of 5 pcs | | 22633 |
| N° 4 | cellpack of 5 pcs | | 22634 |
| N° 1-4 | assort. | box of 24 pcs, 4 reamers | 22595 |
| N° 1-4 | assort. | box of 12 pcs | 22600 |
| N° 1 | box of 12 pcs | | 22596 |
| N° 2 | box of 12 pcs | | 22597 |
| N° 3 | box of 12 pcs | | 22598 |
| N° 4 | box of 12 pcs | | 22599 |

| PD™ Steel abrasive strips | | | REF 0.15 mm 1-sided | REF 0.25 mm 2-sided |
|---------------------------|--------|----------------|---------------------------|---------------------------|
| 2 mm | medium | pack of 12 pcs | 46002 | |
| 3 mm | medium | pack of 12 pcs | 46012 | |
| 4 mm | fine | pack of 12 pcs | 46021 | 46121 |
| 4 mm | medium | pack of 12 pcs | 46022 | 46122 |
| 4 mm | coarse | pack of 12 pcs | 46023 | 46123 |
| 6 mm | fine | pack of 12 pcs | 46031 | |
| 6 mm | medium | pack of 12 pcs | 46032 | 46132 |
| 6 mm | coarse | pack of 12 pcs | 46033 | |
| 8 mm | medium | pack of 12 pcs | 46042 | |

| PD™ PELLA transparent crown forms | | REF |
|-----------------------------------|---|--------------|
| Complete box | 78 crowns (39 pairs) | 33410 |
| Assort. | box of 18 pairs, R&L 1, 11, 21, 31, 41, 51, 101, 111, 121, 131, 2, 12, 22, 32, 3, 13, 23, 5 | 33405 |
| N° 1 | box of 6 pairs | 33421 |
| N° 2 | box of 6 pairs | 33431 |
| N° 3 | box of 6 pairs | 33441 |
| N° 4 | box of 6 pairs | 33451 |
| N° 5 | box of 6 pairs | 33461 |
| N° 6 | box of 6 pairs | 33471 |
| N° 7 | box of 6 pairs | 33481 |
| N° 11 | box of 6 pairs | 33491 |
| N° 12 | box of 6 pairs | 33501 |
| N° 13 | box of 6 pairs | 33511 |
| N° 14 | box of 6 pairs | 33521 |
| N° 15 | box of 6 pairs | 33531 |
| N° 16 | box of 6 pairs | 33541 |
| N° 17 | box of 6 pairs | 33551 |
| N° 21 | box of 6 pairs | 33561 |
| N° 22 | box of 6 pairs | 33571 |
| N° 23 | box of 6 pairs | 33581 |
| N° 24 | box of 6 pairs | 33591 |
| N° 31 | box of 6 pairs | 33601 |
| N° 32 | box of 6 pairs | 33611 |
| N° 33 | box of 6 pairs | 33621 |
| N° 41 | box of 6 pairs | 33631 |
| N° 42 | box of 6 pairs | 33641 |
| N° 43 | box of 6 pairs | 33651 |
| N° 51 | box of 6 pairs | 33661 |
| N° 101 | box of 6 pairs | 33671 |
| N° 111 | box of 6 pairs | 33681 |
| N° 121 | box of 6 pairs | 33691 |

| PD™ PELLA transparent crown forms | | REF |
|-----------------------------------|----------------|--------------|
| N° 131 | box of 6 pairs | 33701 |
| N° 141 | box of 6 pairs | 33711 |
| N° 151 | box of 6 pairs | 33721 |
| N° 161 | box of 6 pairs | 33731 |
| N° 171 | box of 6 pairs | 33741 |
| N° 201 | box of 6 pairs | 33751 |
| N° 202 | box of 6 pairs | 33761 |
| N° 203 | box of 6 pairs | 33771 |
| N° 211 | box of 6 pairs | 33781 |
| N° 212 | box of 6 pairs | 33791 |
| N° 213 | box of 6 pairs | 33811 |

| PD™ Articulating paper | | | REF THIN | REF THICK |
|------------------------|----------|----------------------------|--------------|--------------|
| Strips | blue-red | 12 books of 12 sheets each | 36405 | 36406 |
| Strips | blue | 12 books of 12 sheets each | 36401 | 36402 |
| Strips | red | 12 books of 12 sheets each | 36403 | 36404 |
| Half horseshoe | blue-red | 12 books of 12 sheets each | 36425 | 36426 |
| Half horseshoe | blue | 12 books of 12 sheets each | 36421 | 36422 |
| Half horseshoe | red | 12 books of 12 sheets each | 36423 | 36424 |
| Horseshoe | blue-red | 6 books of 12 sheets each | 36415 | 36416 |
| Horseshoe | blue | 6 books of 12 sheets each | 36411 | 36412 |
| Horseshoe | red | 6 books of 12 sheets each | 36413 | 36414 |

| PD™ Jiffy tubes | | REF |
|-----------------|----------------------------------|--------------|
| Assort. | box of 18 straight and 12 curved | 10682 |
| Straight | box of 36 pieces | 10681 |
| Curved | box of 24 pieces | 10680 |

SPECIAL SUPPLIES CLINIC

PD™ ProphyRing Set REF

| | |
|--|-------|
| Box of 100 disposable cups and 2 ProphyRings | 47080 |
|--|-------|

P-One REF

| | |
|---|-------|
| Box includes 1 dispenser, 1 large grid and 1 small grid | 33280 |
|---|-------|

PD™ Dispenser for cotton pellets REF REF

| | REF smallgrid | REF largegrid |
|-----------------------------------|---------------|---------------|
| Model A standard refill design | 33210 | 33211 |
| Model B pop-up spring bulk design | 33230 | 33240 |

PD™ Dappen dish REF

| | | |
|-------------|---------------|-------|
| Yellow | 1 dappen dish | 11781 |
| Blue | 1 dappen dish | 11782 |
| Green | 1 dappen dish | 11783 |
| Transparent | 1 dappen dish | 11780 |

PD™ Saliva ejectors REF

| | | | |
|--------|--------|-----------------|-------|
| Type A | 155 mm | pack of 100 pcs | 65311 |
| Type B | 125 mm | pack of 100 pcs | 65312 |

PD™ EZ Dropper

| | | REF |
|-------------|---|-------|
| Small | 3 bottles, 6 dappen dishes, 1 tray | 11790 |
| Regular | 4 bottles, 8 dappen dishes, 1 tray | 11791 |
| Premium A | 3 bottles, 6 dappen dishes, 1 dispenser model A, small grid, 1 tray | 11786 |
| Premium A | 3 bottles, 6 dappen dishes, 1 dispenser model A, large grid, 1 tray | 11787 |
| Premium B | 3 bottles, 6 dappen dishes, 1 dispenser model B, small grid, 1 tray | 11788 |
| Premium B | 3 bottles, 6 dappen dishes, 1 dispenser model B, large grid, 1 tray | 11789 |
| Drop bottle | 15 ml, 1 pce, pink | 11760 |
| Drop bottle | 15 ml, 1 pce, yellow | 11761 |
| Drop bottle | 15 ml, 1 pce, blue | 11762 |
| Drop bottle | 15 ml, 1 pce, green | 11763 |

PD™ Cotton pellets REF

| | | REF |
|------|--|-------|
| N° 1 | Ø 8 mm refill (model A), 6 x 3 g, (~630 pcs) | 33251 |
| N° 2 | Ø 6 mm refill (model A), 6 x 3 g, (~1,470 pcs) | 33252 |
| N° 3 | Ø 4.3 mm refill (model A), 6 x 3 g, (~2,220 pcs) | 33253 |
| N° 4 | Ø 3.6 mm refill (model A), 6 x 3 g, (~4,050 pcs) | 33254 |
| N° 5 | Ø 3.2 mm refill (model A), 6 x 3 g, (~6,210 pcs) | 33255 |
| N° 1 | Ø 8 mm bulk (model B), 32 g | 33261 |
| N° 2 | Ø 6 mm bulk (model B), 32 g | 33262 |
| N° 3 | Ø 4.3 mm bulk (model B), 18 g | 33263 |
| N° 4 | Ø 3.6 mm bulk (model B), 18 g | 33264 |
| N° 5 | Ø 3.2 mm bulk (model B), 18 g | 33265 |

PD™ Napkin holders REF

| | | |
|------------------------------|-------|-------|
| Transparent | 40 cm | 65801 |
| Red | 40 cm | 65802 |
| Blue | 40 cm | 65803 |
| White | 40 cm | 65804 |
| Green | 40 cm | 65805 |
| Yellow | 40 cm | 65806 |
| Violet | 40 cm | 65807 |
| Ball chain chrome regular | 40 cm | 65810 |
| Ball chain chrome adjustable | | 65811 |

PD™ Bur cleaning brush REF

| | |
|------------------------|-------|
| Metal slide brass wire | 69302 |
|------------------------|-------|

| PD™ Marking rings | | REF Ø 6.8 mm | REF Ø 8 mm | REF Ø 10 mm |
|-------------------|---------------|-----------------|---------------|----------------|
| Assort. | box of 48 pcs | 21351 | 21371 | 21391 |
| White | box of 48 pcs | 21341 | 21361 | 21381 |
| Yellow | box of 48 pcs | 21342 | 21362 | 21382 |
| Red | box of 48 pcs | 21343 | 21363 | 21383 |
| Blue | box of 48 pcs | 21344 | 21364 | 21384 |
| Green | box of 48 pcs | 21345 | 21365 | 21385 |
| Black | box of 48 pcs | 21346 | 21366 | 21386 |
| Grey | box of 48 pcs | 21347 | 21367 | 21387 |
| Purple | box of 48 pcs | 21348 | 21368 | 21388 |
| Brown | box of 48 pcs | 21349 | 21369 | 21389 |
| Orange | box of 48 pcs | 21350 | 21370 | 21390 |

SPECIAL SUPPLIES LAB

| PD™ Universal solder . Silver base | | REF |
|------------------------------------|-------------------|-------|
| White | pack of 12 sticks | 10001 |

| PD™ Micro-torch | | REF |
|-----------------|--|-------|
| Complete device | | 10036 |

INDEX

| | | | | | |
|-------------------------------------|----|---|----|------------------------------------|----|
| DRILL-A-POST | 59 | PD™ Gutta-percha points . 2 % | 48 | PD™ RubberCup . Soft/Red . Webbed | 18 |
| EssenSeal® | 47 | PD™ Gutta-percha points . 4 % and 6 % | 49 | PD™ Saliva ejectors | 75 |
| FIBRAPOST LUX | 60 | PD™ Gutta-percha points . Conventional | 49 | PD™ Steel abrasive strips | 67 |
| FIBRAPOST PLUS | 61 | PD™ Gutta-percha stopping . Sticks | 54 | PD™ Tofflemire matrix bands | 57 |
| Irriflex® | 28 | PD™ Hold | 54 | PD™ Tofflemire matrix retainers | 57 |
| MAP Accessories | 44 | PD™ Interdental Wooden wedges | 55 | PD™ Transparent rings | 58 |
| MAP Intro Kit | 41 | PD™ Ionover Filling | 64 | PD™ Transparent strips | 58 |
| MAP Needles | 40 | PD™ Ionover Luting | 64 | PD™ Universal solder . Silver base | 78 |
| MAP One | 39 | PD™ Jiffy tubes | 66 | PD™ Vaseline | 56 |
| MAP Surgical Kit | 43 | PD™ Marking rings | 77 | PD™ Zinc Oxide . >99 % pure | 37 |
| MAP Universal Kit | 42 | PD™ Micro-torch | 79 | P-One | 71 |
| OPACAL | 35 | PD™ MTA White | 46 | Pro-Pul-Pan | 36 |
| PD™ Acrylic points . 2 % | 50 | PD™ Napkin holders | 75 | TEMPOTEC | 38 |
| PD™ Anatomical Wooden wedges | 55 | PD™ Paper points . 2 % | 32 | | |
| PD™ Articulating paper | 67 | PD™ Paper points . 4 % and 6 % | 33 | | |
| PD™ Bur cleaning brush | 76 | PD™ Paper points . Conventional | 33 | | |
| PD™ Calcium Hydroxide powder | 62 | PD™ Peeso reamers | 25 | | |
| PD™ Cavity Liner . Light cure | 62 | PD™ PELLA transparent crown forms | 65 | | |
| PD™ C&B Zinc Phosphate cement | 63 | PD™ Plastic pipettes . 0.3 ml | 31 | | |
| PD™ Cotton pellets | 74 | PD™ Plastic pipettes . 1 ml | 30 | | |
| PD™ Dappen dish | 73 | PD™ ProphyCup Fire . Firm/Red . Webbed | 11 | | |
| PD™ Dental floss . Nylon | 21 | PD™ ProphyCup Snow . Soft/White . Webbed | 16 | | |
| PD™ Dental floss . Silk | 21 | PD™ ProphyPaste Pro-F ^{1.23} . Fine . Fluoride | 19 | | |
| PD™ Dental sticks | 21 | PD™ ProphyPaste Pro-F ^{1.23} . Medium . Fluoride | 13 | | |
| PD™ Disclosing tablets | 10 | PD™ ProphyPaste Pro-N ¹⁰⁰ . Fine | 20 | | |
| PD™ Dispenser for cotton pellets | 72 | PD™ ProphyPaste Pro-N ¹⁰⁰ . Medium | 15 | | |
| PD™ EDTA 17 % . Gel | 27 | PD™ ProphyRing Set | 15 | | |
| PD™ EDTA 17 % . Solution | 31 | PD™ Pumice powder | 14 | | |
| PD™ Effervescent mouth wash tablets | 10 | PD™ Resin remover | 26 | | |
| PD™ Eugenate desobturator | 26 | PD™ Root canal barbed broaches | 27 | | |
| PD™ Eugenol . >99 % pure | 37 | PD™ Rosen root canal cement | 51 | | |
| PD™ EZ Dropper | 73 | PD™ RubberCup . Firm/White . Serrated | 12 | | |
| PD™ Gates-Glidden drills | 24 | PD™ RubberCup . Firm/White . Webbed | 12 | | |
| | | PD™ RubberCup . Soft/Red . Serrated | 18 | | |



This image shows a blank sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

This image shows a blank sheet of white paper with horizontal ruling lines. The lines are evenly spaced and extend across the width of the page. There are no margins, text, or other markings on the paper.

Product catalog valid as of 1.9.2019

Products photo credit: Studio BBH Solutions Visuelles, Vevey, Switzerland.

Design and production: Produits Dentaires SA and BBH Solutions Visuelles, Vevey, Switzerland.

© Products, specifications and features can be subject to change without prior notice.



Produits Dentaires SA

Rue des Bosquets 18

1800 Vevey

Switzerland

Phone +41 21 921 26 31

info@pd-dental.com

pd-dental.com